



# epoc

## Blood Analysis System

### Product Specifications

Improve outcomes and workflow, while transforming care delivery across clinical pathways. The epoc® Blood Analysis System gives you Lab quality results here and now.

# epoc<sup>®</sup> Blood Analysis System Product Specifications

System Specifications	
<b>System Description</b>	Point of care blood gas, electrolytes, metabolites and hematocrit analyzer
<b>Sample Types</b>	Whole blood: arterial, venous, mixed venous, capillary
<b>Sample Size</b>	92 µL
<b>Time to Result</b>	Approximately 35 seconds
<b>Calibration</b>	Automatically performed prior to every test
<b>Quality Control</b>	Eurotrol GAS-ISE Metabolite control; Eurotrol Hematocrit control
<b>Integrated Bar-code Scanner</b>	Patient ID and Operator ID; 1D and 2D bar-code format
<b>External Interfaces</b>	HL7
<b>Power Requirements</b>	100–240 VAC; 50–60 Hz (Reader)
<b>Battery</b>	Yes, with typical use up to 70 test cards
<b>Safety</b>	IEC 61010-1, IEC 61010-2-81, IEC 61010-2-101, CSA/UL 601 (Reader)
<b>EMC</b>	IEC60601-1-2 (Reader with AC adapter); IEC 61326-1, IEC 61326-2-6
<b>Operating System</b>	MICROSOFT Windows Mobile 6.5 Classic
<b>Communication</b>	Real Time Wireless, LIS/HIS via data management system, POC Informatics Solutions

Test Card	
<b>Storage</b>	15–30°C
<b>Shelf Life</b>	Up to 5 months
<b>Size</b>	3.39 in (L) × 2.13 in (W) × 0.06 in (H) 86 mm (L) × 54 mm (W) × 1.4 mm (H)

System Dimensions		
	Host	Reader
<b>Length</b>	5.78 in 147 mm	8.46 in 215 mm
<b>Width</b>	3.03 in 77 mm	3.35 in 85 mm
<b>Height</b>	1.06 in 27 mm	2 in 51 mm
<b>Weight</b>	12.5 oz 359 g	<1.1 lb 354 g
<b>Display</b>	3.5 in LCD	—

Environmental Requirements		
<b>Operating Temperature</b>	0°C–50°C (Host)	15°C–30°C (Reader)
<b>Humidity</b>	Up to 95% relative humidity, non-condensing	
<b>Barometric Pressure</b>	400–825 mmHg (53.33-110 kPa)	

epoc and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Measured Parameters		
Parameter	Unit of Measure	Measurement Range
pH	pH units	6.5–8.0
pCO <sub>2</sub>	mmHg kPa	5–250 0.7–33.3
pO <sub>2</sub>	mmHg kPa	5–750 0.7–100
Na <sup>+</sup>	mmol/L mEq/L	85–180
K <sup>+</sup>	mmol/L mEq/L	1.5–12.0
Ca <sup>++</sup>	mmol/L mg/dL mEq/L	0.25–4.00 1.0–16.0 0.5–8.0
Cl <sup>-</sup>	mmol/L mEq/L	65–140
Glu	mmol/L mg/dL g/L	1.1–38.5 20–700 0.20–7.00
Lac	mmol/L mg/dL g/L	0.30–20.00 2.7–180.2 0.03–1.80
Crea	mg/dL µmol/L	0.30–15.00 27–1326
Hct	% PCV L/L	10–75 0.10–0.75

Calculated Parameters		
Parameter	Unit of Measure	Measurement Range
cHgb	g/dL mmol/L g/L	3.3–25 2.0–15.5 33–250
cHCO <sub>3</sub> <sup>-</sup>	mmol/L mEq/L	1–85
cTCO <sub>2</sub>	mmol/L mEq/L	1–85
BE(ecf)	mmol/L mEq/L	-30–+30
BE(b)	mmol/L mEq/L	-30–+30
cSO <sub>2</sub>	%	0–100
eGFR	mL/min/ 1.73m <sup>2</sup>	2–60 or >60*
eGFR-a	mL/min/ 1.73m <sup>2</sup>	2–60 or >60*
AGap	mmol/L mEq/L	-14–+95
AGapK	mmol/L mEq/L	-10–+99
A	mmHg kPa	5-800 0.67-106.64
A-a	mmHg kPa	1-800 0.13-106.64
a/A	% fraction	0-100 0-1

Institutions should establish and set their own normal range values.

“The smart card technology of epoc reduces cost and maximizes efficiency in the hospital setting.”<sup>†</sup>

Agarwal, S et al. Evaluation of the analytical performance of the modified enterprise point-of-care blood gas and electrolyte analyzer in a pediatric hospital. Point of Care 13(4): Dec 2014.

\*Values >60 will be reported as >60 mL/min/1.73 m<sup>2</sup>.

†The outcomes obtained by the Siemens customer described here were realized in the customer’s unique setting. Since there is no typical laboratory, and many variables exist, there can be no guarantee that others will achieve the same results.

**Siemens Healthineers Headquarters**  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84 0  
siemens.com/healthineers

**Local Contact Information**  
Siemens Healthcare Diagnostics Inc.  
Point of Care Diagnostics  
2 Edgewater Drive  
Norwood, MA 02062-4637  
USA  
Phone: +1 781-269-3000  
siemens.com/healthineers